



Dongguan Xintai Instrument Co.,Ltd.

PROFESSIONAL INSTRUMENT MANUFACTURER

## Body Infrared Thermometer HT-820D

- ◆ Non-contact measuring temperature accurately.
- ◆ It can reduce the risk of virus spread effectively.
- ◆ reading temperature display quickly.response time is only 0.5s.
- ◆ Accuracy is  $\pm 0.3^{\circ}\text{C}$  ( $\pm 0.6^{\circ}\text{F}$ )
- ◆ 34 measuring datas can be saved.



Normal



Slight



Serious



### Specification

Storage temperature	-10~60°C(14-140°F)
Usage temperature	10~40°C(50~104°F)
Usage battery	3V
Weight	122g
Size	145 x 80 x 40mm

### Technical Parameters

Temperature measuring range	32~43°C(90-109°F)
temperature measuring accuracy	$\pm 0.3^{\circ}\text{C}$ ( $\pm 0.6^{\circ}\text{F}$ )
Surface Temperature measuring range	0~100°C(32-212°F)
Surface temperature measuring accuracy	$\pm 1^{\circ}\text{C}$ ( $\pm 1.8^{\circ}\text{F}$ )
Minimum readings	$\pm 0.1^{\circ}\text{C}$ ( $\pm 0.1^{\circ}\text{F}$ )
Response time	0.5s
Measuring distance	5~ 8cm
Auto power off	7s

Different temperature range different backlight, distinguished quickly.



MADE IN CHINA



### Distributor:

**AWR Smith Process Instrumentation cc**

PO Box 449 Pinegowrie 2123

Tel: 011 782 7279 / 888 7144 Fax: 086 648 7664

e-mail: [sales@instruments.co.za](mailto:sales@instruments.co.za) website: [www.instruments.co.za](http://www.instruments.co.za)